

# Research Updates

May 17, 2012  
Alex Robinson

# Topics

- Data Reader
- Coherent Scattering Sampling
- Incoherent Scattering Sampling
- Adjoint Incoherent Scattering Sampling
- Code Overview
- Sample Problems


# Data Reader

- An EPDL data reader has been written in C++
- While reading in data, it does necessary manipulations to data so that the data can be used with sampling functions
- This data reader also stores the data arrays for use during particle tracking
- Called the Element Class


# Coherent Scattering Sampling

- Persliden's method has been implemented to efficiently sample coherent scattering
- This method requires numerical integration of the atomic form factor data
- Due to my choice in numerical integration, there is an issue that needs to be worked out


# Coherent Scattering Sampling Efficiency




# Coherent Scattering Sampling Testing (0.1 MeV, H)



# Coherent Scattering Sampling Testing (0.01 MeV, H)




# Coherent Scattering Sampling Testing (0.001 MeV, H)






# Coherent Scattering Sampling Testing (0.0001 MeV, H)




# Incoherent Scattering Sampling

- Kahn's and Koblinger's methods have been implemented to efficiently sample incoherent scattering
- Kahn's method uses a rejection scheme
- Koblinger's method uses a probability mixing scheme
- No numerical integration of data is required by these methods
- No discrepancies have been observed during testing


# Incoherent Scattering Sampling Efficiency




# Incoherent Scattering Sampling Testing (5 MeV, H)




# Incoherent Scattering Sampling Testing (1.0 MeV, H)



# Incoherent Scattering Sampling Testing (0.5 MeV, H)



# Incoherent Scattering Sampling Testing (0.05 MeV)




# Adjoint Incoherent Scattering Sampling


- Sampling methods do not exist in the literature
- I developed a method that relies on both probability mixing and rejection sampling
- High efficiencies are observed with this method
- Good agreement with physical data has been observed




# Adjoint Incoherent Scattering Sampling Efficiency




# Adjoint Incoherent Scattering Sampling Testing (5 MeV, H)




# Adjoint Incoherent Scattering Sampling Testing (1 MeV, H)



# Adjoint Incoherent Scattering Sampling Testing (0.25 MeV, H)



# Adjoint Incoherent Scattering Sampling Testing (0.01 MeV, H)



## Treatment Planning Software

```
graph TD; A[Treatment Planning Software] --> B[Transport Module]; A --> C[Optimization Module]; A --> D[Image Data Read Module]; A --> E[User input];
```

### Transport Module

- Adjoint Monte Carlo
- Forward Monte Carlo
- Adjoint Physics
- Forward Physics
- Combinatorial Geometry
- DAG Geometry

### Optimization Module


- Use adjoint information
- Use Greedy Heuristic
- Brachytherapy
  - Permanent sources for prostate cancer
  - Multi-specie source implant
  - High Dose rate for prostate cancer
  - High Dose rate breast cancer
  - Directional sources
- External beam Radiation Therapy
  - Optimization method
  - Adjoint pencil beams


### Image Data Read Module

- MRI
- Ultra sound
- PET


### User input

- Dose constraints
- Prescribed dose
- Source (seed) type
- Multi-source type










Not coded

Coded but not tested


Coded and tested



Not coded

Coded but not tested

Coded and tested



# Initial MCNP5 Benchmarking

## 10cm Spherical Volume, Point Source, Void

	Surface Current	Surface Flux	Cell Flux
MCNP5	$1.00 \pm 0.0$	$7.95775\text{e-}4 \pm 0.0$	$2.38732\text{e-}3 \pm 0.0$
FAPMC	$7.96\text{e-}4 \pm 0.0$	$7.96\text{e-}4 \pm 0.0$	$2.387\text{e-}3 \pm 0.0$







## 10cm Spherical Volume, 5cm Spherical Volume Source, Void

	Surface Current	Surface Flux	Cell Flux
MCNP5	$1.00 \pm 0.0$	$8.40542\text{e-}4 \pm 0.0$	$2.26389\text{e-}3 \pm 2\text{e-}4$
FAPMC	$7.96\text{e-}4 \pm 0.0$	$8.41\text{e-}4 \pm 3.7\text{e-}5$	$2.263\text{e-}3 \pm 2.39\text{e-}4$

## 10cm Spherical Volume, 1.0 MeV Point Source, Water

	Surface Current	Surface Flux	Cell Flux
MCNP5	$9.93\text{e-}1 \pm 1\text{e-}4$	$8.96204\text{e-}4 \pm 4\text{e-}4$	$2.90082\text{e-}3 \pm 4\text{e-}4$
FAPMC	$7.95\text{e-}4 \pm 2.9\text{e-}5$	$9.04\text{e-}4 \pm 4.14\text{e-}4$	$2.9626\text{e-}3 \pm 4.49\text{e-}4$

# Visualization



# Visualization

